



*What does this new
suction-irrigation handle have
that others don't?*

Multifunction Suction-Irrigation System

Economic • Efficient • Ergonomic • Ecological

The five great advantages compared with others...

economic



- ☐ autoclavable
- ☐ durable
- ☐ low long-term costs
- ☐ universal application

efficient




- ☐ interfaces for various attachments
- ☐ large lumen for efficient irrigation and suction
- ☐ can also be used with pressurised irrigation systems

ergonomic



- ☐ user-friendly operation
- ☐ accurate flow control
- ☐ simple handling
- ☐ straightforward cleaning


ecological



- ☐ reusable instrument
- ☐ no waste disposal problems
- ☐ low irrigation fluid consumption with special valve construction
- ☐ suitable for attachment of reusable tubes

multiple functions

- ☐ simultaneous irrigation / suction while using other instruments
- ☐ irrigation / suction with HF connection
- ☐ pressure irrigation






For pressurised irrigation:

"Laparo Pump" 2215

220 - 230 V~, 50 / 60 Hz
incl.

holder for irrigation stand (2215.991),
reusable tube set (8170.223), power supply
cable, 3 m long (2440.03), vacuum tube
from pump to secretion trap (8170.401),
hygiene filter for aspiration
(4171.121)2215.001

Multifunction Suction-Irrigation System



Suction-irrigation handle8385.901
Sealing cap8385.902
Suction-irrigation tube set, reusable,
 with piercing thorn8170.9521
Tube set, reusable,
 special for laparo pump 22158170.2215
Tube set, PVC, with 2 piercing thorns
 disposable, (pack of 10),
 packed singly, sterile4170.225
Suction-irrigation tube set
 with safe-lock®, (pack of 10),
 packed singly, sterile4170.955
Suction-irrigation tube set,
 with 2 piercing thorns, (pack of 10),
 packed singly, sterile4170.956

Suction-irrigation tube, Ø 3.5 mm,
 310 mm WL8391.911
Suction-irrigation tube, Ø 5 mm,
 310 mm WL8383.911
 450 mm WL8384.911
 310 mm WL,
 curved distal tip (not illustrated) ...8383.912
Pressure irrigation tube, Ø 5 mm,
 310 mm WL8383.913
Suction-irrigation tube, Ø 10 mm, capacity
 5 mm, comprising outer and inner tube,
 225 mm WL8385.942
 310 mm WL8385.951
 required for 8385.942 / 8385.951
Sealing attachment with opening
 and membrane valve8385.903
 Spare membrane valve8920.911

HF spatula probe, Ø 5 mm, monopolar,
 with suction and irrigation channel,
 310 mm WL8383.811
HF hook electrode, Ø 5 mm, monopolar,
 with suction and irrigation channel,
 310 mm WL8383.812
HF button electrode, Ø 5 mm, monopolar,
 with suction and irrigation channel,
 310 mm WL8383.813
Suction-irrigation tube by Bueß,
 insulated, with atraumatic dissection tip
 Ø 10 mm, 330 mm WL8385.944

as an alternative to the suction irrigation handle:
Irrigation attachment
 with two-way tap8383.819

RICHARD WOLF GmbH · 75434 Knittlingen · PF 1164 · Phone +49 7043 - 35-0 · Fax +49 7043 - 3 53 00 · GERMANY · info@richard-wolf.com · www.richard-wolf.com

AUSTRIA:
 RICHARD WOLF Austria
 Ges.m.b.H.
 Wilhelminenstraße 93 a
 A-1160 Wien
 Tel.: 01-405 51 51
 Fax: 01-405 51 51 45
 info@richard-wolf.at

BELGIUM:
 N.V. Endoscopie
 RICHARD WOLF Belgium S.A.
 Industriezone Drongen, Landegemstraat 6
 B-9031 Gent-Drongen
 Tel.: 09.280.81.00
 Fax: 09.282.92.16
 endoscopy@richard-wolf.be

FRANCE:
 RICHARD WOLF France S.A.R.L.
 Rue Daniel Berger
 Z.A.C. La Neuville
 F-51100 Reims
 Tél.: 03.26.87.02.89
 Fax: 03.26.87.60.33
 endoscopes@richardwolf.fr

UK:
 RICHARD WOLF UK Ltd.
 Waterside Way
 Wimbledon
 SW17 0HB
 Tel.: 020 8944 7447
 Fax: 020 8944 1311
 admin@richardwolf.uk.com

USA:
 RICHARD WOLF
 Medical Instruments Corp.
 353 Corporate Woods Parkway
 Vernon Hills, Illinois 60061
 Tel.: 847-913 1113
 Fax: 847-913 1488
 sales&marketing@richardwolfusa.com

WEB (网址) : www.richard-wolf.cc

Printed on paper based on cellulose which has been bleached without the use of chlorine.

Technical specifications subject to change without notice.

sp@stuetzpartner.de

B 669.V.05.GB.1