



关节镜



PANOVIEW HD内窥镜



	4.0mm/175mm	4.0mm/175mm(3 wide)	2.7mm/185mm
0°	8880.541	8885.441	8672.421
30°	8880.543	8885.443	8672.422
45°	8880.544		
70°	8880.545	8885.445	8672.425
90°	8880.437		

AUTOCLAVE
134° C / 273° F

调焦内窥镜

AUTOCLAVE
134° C / 273° F


	f=35mm	f=18mm
30°/170mm	8883.433	8883.473
70°/170mm	8883.435	8883.475

穿刺鞘、锥



快速连接



360度旋转




360°

89121.0260 固定双冲洗通道鞘，直径6.0mm，工作长度130mm，配合4.0mm内窥镜
 89121.1260 旋转双冲洗通道鞘，直径6.0mm，工作长度130mm，配合4.0mm内窥镜
 89121.1263 旋转双冲洗通道鞘，直径6.3mm，工作长度130mm，配合4.0mm内窥镜
 891221255 旋转双冲洗通道鞘，直径5.6mm，工作长度130mm，配合4.0mm内窥镜
 891221241 旋转双冲洗通道鞘，工作长度177.5mm，配合2.7mm内窥镜

排水鞘、锥




8302.09通道内径2.3mm，长度80mm



8303.09通道内径3.7mm，长度73mm

探钩



8863.80直径4.4mm，长度228mm


8864.80直径5.5mm，长度228mm




8399.901钩长7mm

8399.95钩长4.5mm


解剖刀



8436.101直刀




8436.151玫瑰刀



8436.201钩刀



8436.251香蕉刀



8436.301关节刀

锉刀、刮匙、探针



8436.501 直粗锉



8436.511 下弯锉



8436.521 上弯锉



8436.601 刮勺



8436.701 刮匙



8435.901 3mm



8435.902 4.5mm



8435.903 6mm



89160.1009 9mm



89160.1011 11mm




8435.951 15°探针




8435.963 30°探针



8435.953 35°探针



8435.952 45°探针








8435.966 60°探针




8435.956 60°探针



8435.968 80°探针

ARTHROline 咬切钳**咬切钳、带吸引咬切钳、抓钳****剪刀****小型器械****1.9mm小关节镜****2.4mm小关节镜**


ACL前交叉韧带修复



8874.011/012/013/014定位装置




移植模块



综合工作台



8866.951 取腱器



8866.956 ACL槌




8869.991 推动器




8866.101 ACL冲击装置



8866.001/002准直管




空心钻



8869.911 气动手钻



8866.302 ACL管凿装置



打孔钻



8866.932 冲击测试器

动力刨削系统



2304 动力系统



M5/3刨削手柄 89955.0003

带有控制键，最高转速160000转/分，刀头直径2mm-8mm，带吸引通道




M5/0刨削手柄 89955.0000

最高转速160000转/分，刀头直径2mm-8mm，带吸引通道



S1刨削手柄 89951.0001

带有控制键，最高转速3000转/分，刀头直径2mm-5.5mm，带吸引通道



89950.0001

DRILL M1 动力手柄



89950.0400 配用工具直径1.8mm-3.2mm

89950.0500 配用工具直径0.6mm-1.6mm



89950.0300 配用工具直径0mm-6.35mm



89950.0600 配用工具直径0mm-4mm



摆动锯89950.0200



矢状锯 89950.0100

关节灌注泵



2203 关节灌注泵

器械识别

大液晶屏

图形显示

不侵入式压力监测

操作简便

压力监测范围 0-150mmHg

灌注流量范围 0-1.5升/分