







	<p><b>7330.052 3mm 纤维胆道镜</b></p> <p>上弯 130° 下弯 160°                      器械通道 1.5mm                      工作长度 450mm                      包括：                      8734.606 活检钳                      163.903 测漏器                      163.904 气阀                      7321.911 清洁刷                      7305.782 器械箱</p>
 <p>Abwinklung / deflection</p> <p>7305.781 auf Wunsch / on request</p> <p>7305.782 im Lieferumfang enthalten included</p> <p>7305.001 ab / down ← → auf / up</p> <p>7305.006 auf / up ← → ab / down</p>	<p><b>7305.001 5mm 纤维胆道镜</b></p> <p>上弯 210° 下弯 150°                      器械通道 2.5mm                      工作长度 400mm                      包括：                      829.051 活检钳                      163.903 测漏器                      163.904 气阀                      7321.911 清洁刷                      7305.782 器械箱</p>
 <p>Abwinklung / deflection</p> <p>130° auf / up</p> <p>160° ab / down</p>	<p><b>7331.001 3mm 纤维输尿管肾镜</b></p> <p>上弯 130° 下弯 160°                      器械通道 4Fr.                      工作长度 600mm                      包括：                      829.601 活检钳                      163.903 测漏器                      163.904 气阀                      7321.911 清洁刷</p>
 <p>Abwinklung / deflection</p> <p>7305.781 auf Wunsch / on request</p> <p>7305.782 im Lieferumfang enthalten included</p> <p>7305.001 ab / down ← → auf / up</p> <p>7305.006 auf / up ← → ab / down</p>	<p><b>7305.006 5mm 纤维膀胱镜</b></p> <p>上弯 210° 下弯 150°                      器械通道 2.5mm                      工作长度 400mm                      包括：                      7305.782 接合器                      829.061 活检钳</p>

	<p>163. 903 测漏器 163. 904 气阀 7321. 911 清洁刷</p>
 <p>130° 160°</p>	<p><b>7325. 122 2. 5mm 纤维宫腔镜</b></p> <p>上弯 130° 下弯 160° 器械通道 1. 2mm 工作长度 200mm 包括： 7223. 60 活检钳 163. 903 测漏器 163. 904 气阀 7321. 911 清洁刷 7305. 782 器械箱</p>
 <p>210° auf / up 150° ab / down Abwinkelung / deflection</p> <p>7305. 781 auf Wunsch / on request 7305. 782 im Lieferumfang enthalten included</p> <p>7305. 001 ab / down ← → auf / up 7305. 006 auf / up ← → ab / down</p>	<p><b>7305. 011 5mm 纤维输卵管镜</b></p> <p>上弯 210° 下弯 150° 器械通道 2. 5mm 工作长度 400mm 包括： 829. 061 活检钳 163. 903 测漏器 163. 904 气阀 7264. 691 清洁刷 7305. 782 器械箱</p>
	<p><b>7223. 001 3. 5mm 纤维鼻咽喉镜</b></p> <p>上弯 130° 下弯 130° 视野 0° 器械通道 1. 1mm 工作长度 300mm 纤维光导束长 1800mm 包括： 163. 903 测漏器 163. 904 气阀 选配： 7223. 60 活检钳 7223. 65 异物钳 7223. 911 止血电极</p>

	<p><b>7265.001 5mm 纤维支气管镜</b></p> <p>上弯 180° 下弯 130° 视野 115° 器械通道 2.5mm 工作长度 850mm 纤维光导束长 1800mm</p> <p>包括：</p> <ul style="list-style-type: none"><li>7265.851 冲洗接头</li><li>7264.611 活检钳</li><li>163.903 测漏器</li><li>163.904 气阀</li><li>7264.691 清洁刷</li></ul>
	<p>因人体力学及工程制造学获 2004 年医疗器械行业最高奖项—爵士奖的德国著名品牌 WOLF 纤维镜</p>