



Safety first - **The Fluid Manager®**

You operate - we calculate

Straightforward, simple and safe fluid management as a stand-alone solution • Significantly improved efficiency thanks to "Best Flow" automatic instrument detection

The Fluid Manager®

Hystero pump II 2222 + fluid monitoring module 2223



Deficit display

- Compact, clear arrangement of all system components as a stand-alone solution
- Microprocessor-controlled irrigation pump 2222 for hysteroscopy, optional with fluid monitoring module 2223 for intraoperative fluid management
- "Best Flow"
Automatic instrument detection to optimize pressure and flow characteristics, particularly in small-lumen instruments
- In difficult viewing conditions:
"High-Flow" mode (500 ml/min)
"High-Pressure" mode (80-150 mmHg)



Hystero pump II 2222

The Fluid Manager®

Complete system comprising:
Hystero pump II (2222.001) and Fluid monitoring module (2223.001),
100-230 V a.c., 50/60 Hz**2222.023**
with CAN bus interface2222.323

Hystero pump II

incl. tube set, silicone, re-usable (8170.223)
and power cable (2440.03),
100-230 V a.c., 50/60 Hz**2222.001**
with CAN bus interface2222.101

Fluid monitoring module

comprising:
Deficit display, mobile stand with scales and rings to accommodate up to 5 suction containers, device console, splash protection, support beam
for bags2223.001

Further accessories:

Tube set,
silicone, reusable8170.223

Tube set,
with two piercing thorns, sterile,
disposable (pack of 10)4170.223

Tube set,
with safe-lock connector, sterile
disposable (pack of 10)4170.224

Membrane
with strap for pressure chamber
disposable (pack of 10)8170.231

Suction container, 3l,
disposable (pack of 2)2215.971

Suction container, 3l,
autoclavable8170.981

Ring
to support 2215.971
(supplied with 2223.001)64330.284

Fluid collecting apron,
sterile, disposable
(pack of 5)4170.910
also:

Tube set,
Outlet, sterile, disposable
(pack of 10)4170.901

RICHARD WOLF GmbH · 75434 Knittlingen · PF 1164 · Phone +49 70 43 35-0 · Fax +49 70 43 3 53 00 · GERMANY · info@richard-wolf.com · www.richard-wolf.com

AUSTRIA:
RICHARD WOLF Austria
Ges.m.b.H.
Wilhelminenstraße 93 a
A-1160 Wien
Tel.: 01-405 62 19
Fax: 01-405 51 51 45
info@richard-wolf.at

BELGIUM:
N.V. Endoscopie
RICHARD WOLF Belgium S.A.
Industriezone Drongen, Landegemstraat 6
B-9031 Gent-Drongen
Tel.: 09.280.81.00
Fax: 09.282.92.16
endoscopy@richard-wolf.be

FRANCE:
RICHARD WOLF France S.A.R.L.
Rue Daniel Berger
Z.A.C. La Neuville
F-51100 Reims
Tél.: 03.26.87.02.89
Fax: 03.26.87.60.33
endoscopes@richardwolf.fr

UK:
RICHARD WOLF UK Ltd
Waterside Way
Wimbledon
SW17 0HB
Tel.: 020 8944 7447
Fax: 020 8944 1311
admin@richardwolf.uk.com

USA:
RICHARD WOLF
Medical Instruments Corp.
353 Corporate Woods Parkway
Vernon Hills, Illinois 60061
Tel.: 847-913 1113
Fax: 847-913 1488
sales&marketing@richardwolfusa.com

TEL : 22158533/22159298 WEB : www.richard-wolf.cc